



Single-Photon Software Lock-in for PL

Pipelines, Results, and Applications

Summary

1. Built two photon-counting lock-in pipelines:
 - Jakob et al. + Liu et al., yields a photon count.
 - Jakob et al. + Braun & Libchaber, yields modulation depth M and phase ϕ .
2. Explored different experimental conditions, each time comparing the lock-in and regular counting.
3. Implemented FDLM : extracting τ from $\phi(\omega)$ and/or $M(\omega)$.

Deliverables: CLI tools (reader, lock-ins), simple decision rules for when to use lock-in vs counting, and a (preliminary) FDLM workflow.



Quick Background

Classical:

- Inputs: Signal $V_s(t)$, Reference $V_r(t)$
- Mixing (multiplication) -> low pass filtering (integration)
- Outputs: R , θ

Photon counting:

- Inputs: $\{t_p\}$, $\{t_m\}$ taken in $[0, T]$
- Lock-in routine
- Outputs: Photon count ($\leq \text{card}(\{t_p\})$), mod. depth, θ




Physical lock-in inner workings

Hardware & Raw Data

Setup: Laser → Chopper (out: TTL) → Optics → SPAD (out: NIM) → PicoHarp-330 (T2).

Raw data (.PTU) (look inside, bits) -> decode -> .txt
In practice: .PTU -> lock-in (.txt files are huge)




```
75 -----
76 PicoHarp 300 T2 data
77
78 record# chan  nsync truetime/ps
79 0 MAR  2 3784898
80 1 CHN 1 27790538 111162152
81 2 MAR  2 28787954
82 3 CHN 1 42820706 171282824
83 4 CHN 1 42844464 171377856
84 5 CHN 1 43200575 172802300
85 6 CHN 1 48522661 194090644
86 7 CHN 1 48831309 195325236
87 8 CHN 1 49007922 196031688
88 9 MAR  2 53797442
89 10 CHN 1 67651075 270604300
```

Algorithms


Liu et al. (reference-weighted counting) and Braun & Libchaber

Algorithm A — Liu (reference-weighted counting)



- Compute instantaneous phase for each photon from reference markers.
- Accumulate $I = \sum \text{square}(\phi)$, $Q = \sum \text{square}(\phi + 90^\circ)$ over windows
- Signal Strength (or Photon Count) = $\text{norm}(I + iQ)$, Phase = $\text{atan2}(Q, I)$
- Jakob Augmentation: dynamic references, better per-photon phase calculation.
- Main use cases: Multiple sources and responses, daylight operations.

Corrected Signal Modulation Phase Histogram



Algorithm A — Braun & Libchaber

Algorithm 1 Photon-Counting Lock-in Detection

```

1: function PHOTONCOUNTINGLOCKIN

2:   // Step 1: Record photon counts in time bins
3:    $N_{raw} \leftarrow$  empty array
4:   for each time bin  $t_n$  with interval  $1/f_{sample}$  do
5:      $N_{raw}[n] \leftarrow$  number of photons detected in bin
6:   end for

7:   // Step 2: Create functional basis
8:    $R_{raw\_avg} \leftarrow \frac{1}{T} \sum_{n=0}^{T-1} R_{raw}(t_n)$  ▷ Time average
9:   Find zero crossings  $t_z$  with positive slope
10:   $f_{ref}, \phi_{ref} \leftarrow \text{LinearRegression}(t_z)$  ▷ Extract frequency and phase
11:   $R_X(t_n) \leftarrow \sin(2\pi f_{ref} t_n + \phi_{ref})$  ▷ Sine reference
12:   $R_Y(t_n) \leftarrow \cos(2\pi f_{ref} t_n + \phi_{ref})$  ▷ Cosine reference

13:  // Step 3: Project counts onto basis
14:   $I \leftarrow \frac{1}{T} \sum_{n=0}^{T-1} N_{raw}(t_n)$  ▷ Average count rate
15:   $N(t_n) \leftarrow N_{raw}(t_n) - I$  ▷ Center around zero
16:   $I_X \leftarrow \frac{\sum_n N(t_n) \cdot R_X(t_n)}{\sqrt{\sum_n R_X(t_n)^2}}$  ▷ In-phase component
17:   $I_Y \leftarrow \frac{\sum_n N(t_n) \cdot R_Y(t_n)}{\sqrt{\sum_n R_Y(t_n)^2}}$  ▷ Quadrature component

18:  // Step 4: Calculate amplitude and phase
19:   $A \leftarrow \frac{\sqrt{I_X^2 + I_Y^2}}{I}$  ▷ Relative amplitude
20:   $\phi \leftarrow \arctan \frac{I_Y}{I_X}$  ▷ Signal phase
21:  return  $A, \phi, I$ 

22: end function

```





Fig. 1. Implementation of the photon-counting lock-in. Projections of time-binned counts N_{raw} to sinusoidal reference R_X and cosinusoidal reference R_Y are time averaged to yield the small I_X and I_Y signals of the lock-in.



Jakob et al. Augmentation: no need to bin, we resolve each individual photon arrival.

Baselines & Statistics

OFF estimates background; subtract from ON.


SNR = mean/std







SNR Results

SNR — Lock-in vs ON/OFF





In photon counting without locking, we can measure only photon counts I at the noise limit [Eq. (4)]. When we use the photon-counting lock-in, we can also detect the small-signal amplitude and small-signal phase at the same fundamental detection limit [Eqs. (5)] as shown in Fig. 3.

Some Applications

Multiplexing, Lifetimes from Phase and Modulation







FIGURE 1. Schematic of the fNIRS-DOT instrument.



FDLM Theory (single-exponential)




As modulation frequency is increased:



- Modulation decreases from 1 to 0
- Phase increases from 0, $\pi/2$

Idea: Extract lifetime from dynamics $M(\omega)$ and $\phi(\omega)$. Models are especially simple if we assume single exponential decay:


- $M(\omega) = 1 / \sqrt{1 + (\omega\tau)^2}$
- $\tan \phi(\omega) = \omega \tau$



Use phase and/or modulation to fit τ ; amplitude a bit more straightforward.



Fit of $1 / \sqrt{1 + (\tau\omega)^2}$



FDLM Pipeline

- 1) Frequency sweep: collect TTTR T2 at $f_1 \dots f_n$.
- 2) Braun(f) \rightarrow (M , ϕ).
- 3) Fit τ from $M(\omega)$.


```
● └─[$] <git:(main*)> python3 braun_lockin_original.py data/GeSb/10mW_GeSb_10s_10KHz_50pctduty cycle_000.ptu  
File: data/GeSb/10mW_GeSb_10s_10KHz_50pctduty cycle_000.ptu  
RESULT – Braun photon-counting lock-in:  
  f_ref           : 9999.860445 Hz  
  f_sample        : 499993.022250 Hz   (M = 50 bins/period)  
  bins / photons  : 4999797 / 563984  
  I_dc (mean/bin) : 0.112801 counts/bin  
  A (mod depth)   : 1.21884303 (dimensionless)  
  phi             : 3.118130083 rad   (178.655694 deg)  
  bins_per_period : 50  
  trimmed periods : 1 at start and end
```

Lessons Learned, Recs, If I had more time...


- Lock-in is not a free lunch: with stable background, ON–OFF sets a high bar in terms of SNR.
- Lock-in enables measurements beyond frequency counting: modulation, phase, multiple responses, daylight measurements (!)
- Use lock-in when we know of periodic noise, nearby in the band.
- When we need to enable measurements stated above.
- Otherwise, not worth the overhead. Introduces extra time and complexity (sources of error).
- Multi-exponential FDL fits, and further improvements to the models (regularization, error propagation, etc.)
- Cross-validation against TCSPC decay fits
- Dual-modulation experiment ?!

- Jakob, L. A., Deacon, W. M., Hicks, O., Manyakin, I., Ojambati, O. S., Traxler, M., & Baumberg, J. J. (2021). *Single-photon multiclock lock-in detection by picosecond timestamping*. **Optica**, 8(12), 1646–1653.
<https://doi.org/10.1364/OPTICA.441487>
- Braun, D., & Libchaber, A. (2002). *Computer-based photon-counting lock-in for phase detection at the shot-noise limit*. **Optics Letters**, 27(16), 1418–1420.
<https://doi.org/10.1364/OL.27.001418>
- Liu, D., Wang, B., Pan, T., Li, J., Qin, Z., Zhang, L., Zhou, Z., & Gao, F. (2015). *Toward quantitative near infrared brain functional imaging: Lock-in photon counting instrumentation combined with tomographic reconstruction*. **IEEE Transactions on Biomedical Engineering**, 62(12), 2941–2951.
<https://doi.org/10.1109/TBME.2015.2450217>
- PicoQuant GmbH. (2013). *PicoHarp 300: Picosecond histogram accumulating real-time processor — User's manual and technical data* (Software version 2.3). Berlin, Germany: PicoQuant GmbH.
- Lakowicz, J. R. (2006). *Principles of Fluorescence Spectroscopy* (3rd ed.). Springer. (Frequency-Domain Lifetime Measurements, Ch. 5).

Reconstructed Laser Modulation (Run 15; 750 Hz Modulation). Bin Size = $5.33\text{e-}06\text{s}$



Typical Reconstructed ON/OFF Run (PicoHarp 300 data)



Lock-In Outputs at Different Laser Powers, High and Low Background.

Signal vs Optical Density - Linearity Test

